

Intranasal Self-Administration SOP

Developed by the [Siciliano Lab](#)

Department of Pharmacology, Vanderbilt University

Contents

Materials	2
General Experimental Guidelines	2
Rolling Cap Session Setup Instructions	4

Materials

See <https://github.com/Siciliano-Lab/INSA> for a complete list of parts. Note that any standard head-fixed operant setup can be used in place of the described setup; the only parts that are not found in a standard setup are the high resolution (sub-microliter) syringe pump, miniature micromanipulator, and blunt needle.

General Experimental Guidelines

1. Preparation
 - a. Transport animals to the behavior room and allow them to habituate to the room and ambient noise for ≥ 30 minutes.
 - b. In the behavior room:
 - i. Turn on all the necessary equipment including power cabinet, MED-PC V software, lever, pump, and red light inside the chamber.
 - ii. Prepare the video recording system.
2. Prime the infusion line
 - a. Prepare the cocaine solution (concentration as per experimental design).
 - b. Load the syringe and attach it to the microinfusion pump.
 - c. Fill the tubing with cocaine solution until fluid nearly exists the blunt needle.
 - d. Manually advance the syringe by turning the pump knob to ensure consistent droplet formation at the needle tip.
 - i. The volume of cocaine depends on tubing length, diameter, and syringe size.
3. System check
 - a. Open MED-PC V and load the appropriate behavioral program.
 - b. Confirm input detection by checking that the lever presses register on the live input display.
 - i. If no input is detected, check wiring and device connections.
 - c. Enter subject ID and session details.
 - d. Ensure that each triggered delivery reliably results in fluid protruding from the needle
 - e. *If running the rolling cap MED-PC code, see [below](#).*
4. Box set-up
 - a. Wipe down all the apparatus components (tunnel, platform, etc) with 0.03% acetic acid (or preferred cleaning solution).
 - b. Secure tunnel to platform.
5. Head-fixation
 - a. Gently lift the mouse by the tail and allow them to maneuver into the tunnel. Maintain a grip on their tail.
 - b. Gently hold the mouse by one side of the implanted headbar and inset the otherside into the holder and tighten the screw.
 - i. You can now let go of the tail.
 - c. Attache and secure the opposite side of the headbar.
 - d. Ensure:
 - i. The mouse is centered on the platform.
 - ii. The forelimbs can comfortably reach the lever.
 - iii. The body is not too low or overextended – adjust platform height if needed.
6. Position the infusion needle:
 - a. Using the 3-axis micromanipulator, align the blunt needle with the mouse's nostril (**Figure 1**).
 - i. Position the needle tip ~ 1 mm away from the nostril opening.
7. Run the behavioral session

- a. Start the video recording.
 - b. Start the MED-PC behavioral program.
 - i. Again, if running rolling cap program, ensure the macros have been configured.
 - c. Monitor the session until the program terminates or animal displays any signs of distress.
8. Post-session procedure
 - a. Stop the video.
 - b. Swivel the delivery needle away from the mouse.
 - c. Unscrew both ends of the headbar holder and remove the mouse.
 - d. Weigh the animal and return to its home cage.
9. Clean-up
 - a. Wipe down the platform and headfix set up with 0.03% acetic acid (or preferred cleaning solution).
 - b. Remove the syringe from the pump.
 - c. Expel any remaining cocaine into the appropriate biohazard container.
 - d. Flush tubing (including blunt needle) with water followed by air to remove any residual drug solution.
 - e. Turn off all equipment and transfer data.

Rolling Cap Session Setup Instructions

Follow the additional instructions below to execute a self-administration session in which mice receive intranasal cocaine deliveries with a *rolling cap*. To prevent overdose, the rolling cap programs are written such that the reinforcement schedule is temporarily discontinued if the rolling cap (e.g. 40 mg/kg) is reached within the timeframe specified (e.g. 60 minutes). The session continues regardless, and the schedule is re-instated as soon as the specified condition is false. To run these programs, it is necessary to specify the cocaine limit (in mg/kg), the rolling window duration (in minutes), and the cocaine dose that is being delivered (in mg/kg/delivery) through the use of a macro file.

1. Pre-session set up (only performed once per computer)
 - a. Power on the MED-PC cabinet and open MED-PC V.
 - b. In the menu bar navigate to:
Macros > Configure Context Dependent Macros > Add Program, then click OK.
 - i. Select and add the rolling cap program from the list.
 - c. With the rolling cap program highlighted, click Add macros to program, select CocaineDose.MAC, then click OK
 - d. Click OK again to confirm and exit setup
2. Test Session Initialization
 - a. On the main page select the appropriate box and load rolling cap program.
 - b. With the box selected, go to context macros, and choose CocaineDose.mac
 - c. When prompted, enter the following parameters:
 - i. Cocaine Dose (mg/kg/delivery)
 - ii. Window size (minutes)
 - iii. Cocaine Cap (mg/kg)
 - d. Test the program.
3. Follow the general instructions described in the section [above](#)